

## Vala®Protect Basic, Special and Special Eco



Spec.-No.: D 9.0301  
 Department: COH-IDT  
 Date: 2018-04-27

### 1. General Product Description

Vala®Protect disposable protective bed sheets with absorbent tissue and impermeable plastic foil underside, Vala®Protect Special and Special Eco are reinforced with yarn and therefore particularly hard-wearing and tear resistant.

### 2. Application / Indication

Vala®Protect disposable protective bed sheets protect patients, bed linen and examination beds from contamination.

### 3. Presentations

| Art. No. |                                   | package |      |
|----------|-----------------------------------|---------|------|
| 992 226  | Vala®Protect Basic 38 x 38        | P 250   |      |
| 992 227  | Vala®Protect Basic 80 x 140       | P 100   |      |
| 992 228  | Vala®Protect Basic 80 x 175       | P 25 *  | OK 4 |
| 992 229  | Vala®Protect Basic 80 x 210       | P 100   |      |
| 992 238  | Vala®Protect Special 80 x 175     | P 100   |      |
| 992 239  | Vala®Protect Special 80 x 210     | P 100   |      |
| 992 258  | Vala®Protect Special Eco 80 x 175 | P 100   |      |
| 992 259  | Vala®Protect Special Eco 80 x 210 | P 100   |      |

\* - *pharmacy pack*

## Vala®Protect Basic, Special and Special Eco



Spec.-No.: D 9.0301

Department: COH-IDT

Date: 2018-04-27

### 4. Product Characteristics

**Vala®Protect Basic**                      tissue, TCF bleached  
PE-film, green  
adhesive

**Vala®Protect Special**                    tissue, TCF bleached  
PE-film, green  
adhesive  
threads, white

**Vala®Protect Special Eco**                tissue, unbleached  
PE – film, green  
adhesive  
threads, white

### 5. Product Requirements

|                                     | <b>total length<br/>cm</b>              | <b>total width<br/>cm</b>            | <b>weight<br/>g</b>                   |
|-------------------------------------|---|--------------------------------------|---------------------------------------|
| <b>Vala®Protect Basic</b>           | 38 ± 2<br>140 ± 6<br>175 ± 6<br>210 ± 6 | 38 ± 2<br>78 ± 4<br>78 ± 4<br>78 ± 4 | 5 ± 1<br>48 ± 6<br>60 ± 10<br>71 ± 10 |
| <b>Vala®Protect Special</b>         | 175 ± 6<br>210 ± 6                      | 78 ± 4<br>78 ± 4                     | 70 ± 10<br>85 ± 10                    |
| <b>Vala®Protect Special<br/>Eco</b> | 175 ± 6<br>210 ± 6                      | 78 ± 4<br>78 ± 4                     | 75 ± 10<br>90 ± 11                    |

## Vala®Protect Basic, Special and Special Eco



Spec.-No.: D 9.0301  
Department: COH-IDT  
Date: 2018-04-27

6. Labelling

Lot-No. with 9-digit code:  
e.g.:

|      |     |      |               |
|------|-----|------|---------------|
| X    | X   | XX   | XXXXX         |
| year | day | week | internal code |

7. Packaging

Vala®Protect Special, Special Eco and Basic (except for 38 x 38 and 80 x 175 pharmacy pack) packaged with banderole, transport carton acc. DIN, sealed with adhesive tape on europallet.

Vala®Protect Basic 38 x 38 packaged in PE-bag; transport carton acc. DIN, sealed with adhesive tape on europallet.

Vala®Protect Basic 80 x 175 pharmacy pack packaged in PE-bag; transport carton acc. DIN, sealed with adhesive tape on europallet.

April 27<sup>th</sup>, 2018

PAUL HARTMANN AG  
- Product Development COH-IDT -

i.A.

